


SPECIFICATIONS

2,3-BIS(2-METHOXY-4-NITRO-5-SULFOPHENYL)-5-[(PHENYLAMINO)CARBONYL]-2H-TETRAZOLIUM INNER SALT, SODIUM SALT (1:1)



Catalogue Number	2525
CAS Registry Number	111072-31-2
Molecular Formula	C₂₂H₁₆N₇O₁₃S₂Na
Molecular Weight	673.52 g.mol⁻¹
O_{max}²⁸⁵	15,800 ± 1000
¹H- NMR	Consistent with assigned structure
% Water	< 15%
Appearance	Yellow solid
Storage	2-8°C
Use	Reagent for use in cell culture assay systems