



SPECIFICATIONS

cGMP DEHYDROASCORBIC ACID DIMER (DHAA)

Appearance	Off-white to beige colored solid
Solubility	Soluble
Identification by IR	Conforms to reference
Identification by 13C NMR	Conforms to reference
Water Content	Report Only (%)
Functional Assay by UV	NLT 80%
Residue on Ignition	Report Only
Residual Solvents by GC-HSS	
Methanol	Report Only
Acetone	Report Only
Methyl ethyl ketone (MEK)	Report Only
Product Use	Used as a thiol oxidizing agent.
Storage	15 - 30°C under Nitrogen. Protect from water/humidity