


SPECIFICATIONS

DITHIOERYTHRITOL



Catalogue Number	1360
CAS Registry Number	6892-68-8
Molecular Formula	C₄H₁₀O₂S₂
Molecular Weight	154.25 g.mol⁻¹
Melting Point	> 81EC
Titration ex -SH (iodometric)	∃ 99.0%
T.L.C.	One Spot
Appearance	White crystals
5% Solution in Water	Clear and colorless
A₂₈₃ (0.02 M Aqueous Solution)	# 0.040
Storage	2-8EC under Nitrogen