




SPECIFICATIONS

p-NITROPHENYLPHOSPHATE, DISODIUM SALT HEXAHYDRATE



Catalogue Number	2370
CAS Registry Number	4264-83-9
Molecular Formula	C₆H₄NO₆PNa₂·6H₂O
Formula Weight	371.1 g.mol⁻¹
0max³¹²	≧ 9,600
p-Nitrophenol	# 0.005%
% Water	30% ∇ 1.5%
Appearance	White to pale yellow crystalline solid
Percent Free Inorganic Phosphate	# 1.0%
Δ OD (14 Day Stability Test)	# 0.025 AU
Storage	2-8EC