




## SPECIFICATIONS

### p-NITROPHENYLPHOSPHATE, diTRIS SALT



<b>Catalogue Number</b>	<b>2372</b>
<b>CAS Registry Number</b>	<b>68189-42-4</b>
<b>Molecular Formula</b>	<b>C<sub>14</sub>H<sub>28</sub>N<sub>3</sub>O<sub>12</sub>P</b>
<b>Molecular Weight</b>	<b>461.4 g.mol<sup>-1</sup></b>
<b>Omax<sup>312</sup></b>	<b>≥ 9,700</b>
<b>p-Nitrophenol</b>	<b># 0.005%</b>
<b>% Water</b>	<b># 1.0%</b>
<b>Appearance</b>	<b>White to off-white solid</b>
<b>Storage</b>	<b>2-8°C</b>