


SPECIFICATIONS

Page 1 of 1

ISOPROPYL- β -D-THIOGALACTOPYRANOSIDE, NON-ANIMAL ORIGIN

CAT. NO.: 1900



CAS Registry Number	367-93-1
Molecular Formula	C ₉ H ₁₈ O ₅ S
Molecular Weight	238.31 g.mol ⁻¹
Appearance	White crystalline powder
Solubility (200 mM aq)	Clear and colorless
Melting Point	112°C ± 2°C
I.R.	Conforms
Specific Rotation	-31.5 ± 3.0°
Purity	≥ 99.0%
Storage	2-8°C

SPEC1900-1
April 10, 2018