




SPECIFICATIONS

1-(3-DIMETHYLAMINOPROPYL)-3-ETHYLCARBODIIMIDE HYDROCHLORIDE (EDAC-HCl)



Catalogue Number	1339
CAS Registry Number	25952-53-8
Molecular Formula	C₈H₁₇N₃AHCl
Molecular Weight	191.7g.mol⁻¹
Molar Absorptivity	6800-7500
Chloride	18.12 - 18.86% (w/w)
IR	Consistent with assigned structure
Melting Point	108-112°C
Appearance	White, crystalline powder
Storage	-15 ± 5°C under Nitrogen
Use	Peptide Coupling Reagent