


**SPECIFICATIONS****PHOSPHOENOLPYRUVATE TRICYCLOHEXYLAMMONIUM SALT**

<b>Catalogue Number</b>	<b>2555</b>
<b>CAS Registry Number</b>	<b>35556-70-8</b>
<b>Molecular Formula</b>	<b>C<sub>21</sub>H<sub>44</sub>N<sub>3</sub>O<sub>6</sub>P</b>
<b>Molecular Weight</b>	<b>465.6 g.mol<sup>-1</sup></b>
<b>Omax<sup>230</sup></b>	<b>3,000 ± 200</b>
<b>% Water</b>	<b>4% ± 2%</b>
<b>Appearance</b>	<b>White crystalline solid</b>
<b>Storage</b>	<b>2-8°C</b>